


The Home Page has a title bar greeting. A bottom navigation bar is present on all pages. It's destinations are the Home page, the Action Page and the user Profile Page.




The Home Page contains general information regarding the functionality of the application.

The Action Page is a chronological view of all assessments and journal entries. The FABs allow for entering notes, contraction assessments, and fetal assessments respectively.



On the Profile Page, the user enters general demographic information and relevant medical history. The profile can be updated at any time.



The Journal FAB will lead to a pop-up where a timestamped note can be entered. (To be designed)

The Contraction FAB will lead to a page where the onset, pain intensity, and resolution time of a contraction can be entered. The FHR can also be added here. (To be designed)

The FHR FAB will lead to a pop-up where the FHR can be recorded via a USB device attachment. (To be designed)